

AUDIT TRAIL FOR BALANCE

DATE	TIME	USER ID	ACTION	BALANCE CODE NO.	OLD VALUES	NEW VALUES	OLD WEIGHT	NEW WEIGHT	REMARK
16.10.2024	10:19:46	User ID:123,User Name:yogita	Edit	Balance@22	NA	NA	Daily:Standard Weight:10.000,(-)Tolerance:9.500,(+)Tolerance:10.500 Daily:Standard Weight:20.000,(-)Tolerance:18.500,(+)Tolerance:20.500 Removed Data:Standard Weight:100.000,(-)Tolerance:99.000,(+)Tolerance:101.000 Daily:Checked,Eccentricity:Unchecked,Repeatability:Unchecked,Periodic:Unchecked,Uncertainty:Unchecked	Daily:Standard Weight:70.000,(-)Tolerance:69.000,(+)Tolerance:71.000 Daily:Standard Weight:80.000,(-)Tolerance:79.000,(+)Tolerance:81.000 Daily:Standard Weight:90.000,(-)Tolerance:89.000,(+)Tolerance:91.000	

AUDIT TRAIL FOR BALANCE

DATE	TIME	USER ID	ACTION	BALANCE CODE NO.	OLD VALUES	NEW VALUES	OLD WEIGHT	NEW WEIGHT	REMARK
16.10.2024	10:08:55	User ID:123,User Name:yogita	Add	Balance@22	NA	Balance Code No.: NA Balance@22,Make: Sartorius,Model: BC-ED,Serial No.: NA,Unit: gm,Least Count: 0.002,Maximum Capacity: 350.000,Minimum Capacity: 0.002,Maximum Operating Range: 256.000,Minimum Operating Range: 0.003,Periodic Data: Calibration Date:16-10-2024,Duration :5,Reminder:4,Balance Type:Analytical	NA	Daily: Standard Weight:10.000,(-)Tolerance:9.500,(+)Tolerance:10.500 Standard Weight:20.000,(-)Tolerance:18.500,(+)Tolerance:20.500 Standard Weight:30.000,(-)Tolerance:29.000,(+)Tolerance:30.500 Standard Weight:40.000,(-)Tolerance:39.000,(+)Tolerance:41.000 Standard Weight:50.000,(-)Tolerance:40.000,(+)Tolerance:51.000 Standard Weight:60.000,(-)Tolerance:59.000,(+)Tolerance:61.000 Standard Weight:100.000,(-)Tolerance:99.000,(+)Tolerance:101.000 Periodic: Standard Weight:10.000,(-)Tolerance:9.500,(+)Tolerance:10.500 Standard Weight:20.000,(-)Tolerance:18.500,(+)Tolerance:20.500 Standard Weight:30.000,(-)Tolerance:29.000,(+)Tolerance:30.500 Standard Weight:40.000,(-)Tolerance:39.000,(+)Tolerance:41.000 Standard Weight:50.000,(-)Tolerance:40.000,(+)Tolerance:51.000 Standard Weight:60.000,(-)Tolerance:59.000,(+)Tolerance:61.000 Standard Weight:70.000,(-)Tolerance:69.000,(+)Tolerance:71.000 Standard Weight:80.000,(-)Tolerance:79.000,(+)Tolerance:81.000 Standard Weight:90.000,(-)Tolerance:89.000,(+)Tolerance:91.000	test

AUDIT TRAIL FOR BALANCE

DATE	TIME	USER ID	ACTION	BALANCE CODE NO.	OLD VALUES	NEW VALUES	OLD WEIGHT	NEW WEIGHT	REMARK
13.10.2024	10:59:27	User ID:A,User Name:Abhishek	Edit	12345678900987 65432112345	Periodic Data:Duration:999,Reminder:5	Periodic Data:Duration:9,Reminder:9	Periodic:Standard Weight:30.00000,(- Tolerance:29.00000,(+)Tolerance:31. 00000,Periodic:Standard Weight:40.00000,(- Tolerance:39.00000,(+)Tolerance:41. 00000,Periodic:Standard Weight:50.00000,(- Tolerance:49.00000,(+)Tolerance:51. 00000,Periodic:Standard Weight:60.00000,(- Tolerance:59.00000,(+)Tolerance:61. 00000,Periodic:Standard Weight:70.00000,(- Tolerance:69.00000,(+)Tolerance:71. 00000	NA	rt

AUDIT TRAIL FOR BALANCE

DATE	TIME	USER ID	ACTION	BALANCE CODE NO.	OLD VALUES	NEW VALUES	OLD WEIGHT	NEW WEIGHT	REMARK
13.10.2024	10:15:55	User ID:11,User Name:Mangesh	Activate	Balance-121212	NA	NA	NA	NA	9

AUDIT TRAIL FOR BALANCE

DATE	TIME	USER ID	ACTION	BALANCE CODE NO.	OLD VALUES	NEW VALUES	OLD WEIGHT	NEW WEIGHT	REMARK
13.10.2024	10:13:45	User ID:A,User Name:Abhishek	Add	1234567890098765432112345	NA	Balance Code No.: 1234567890098765432112345,Make: METTLER TOLEDO,Model: ML203,Serial No.: NA,Unit: gm,Least Count: 00000.00001,Maximum Capacity: 32000.00000,Minimum Capacity: 00000.00001,Maximum Operating Range: 25000.00000,Minimum Operating Range: 00000.00020,Periodic Data: Calibration Date:28-10-2024,Duration :999,Reminder:5,Balance Type:Analytical	NA	Daily: Standard Weight:0.020,(-)Tolerance:0.019,(+)Tolerance:0.021 Standard Weight:10,(-)Tolerance:9,(+)Tolerance:11 Standard Weight:20,(-)Tolerance:19,(+)Tolerance:21 Standard Weight:30,(-)Tolerance:29,(+)Tolerance:31 Standard Weight:40,(-)Tolerance:39,(+)Tolerance:41 Standard Weight:50,(-)Tolerance:49,(+)Tolerance:51 Standard Weight:60,(-)Tolerance:59,(+)Tolerance:61 Standard Weight:70,(-)Tolerance:69,(+)Tolerance:71 Periodic: Standard Weight:0.020,(-)Tolerance:0.019,(+)Tolerance:0.021 Standard Weight:10,(-)Tolerance:9,(+)Tolerance:11 Standard Weight:20,(-)Tolerance:19,(+)Tolerance:21 Standard Weight:30,(-)Tolerance:29,(+)Tolerance:31 Standard Weight:40,(-)Tolerance:39,(+)Tolerance:41 Standard Weight:50,(-)Tolerance:49,(+)Tolerance:51 Standard Weight:60,(-)Tolerance:59,(+)Tolerance:61 Standard Weight:70,(-)Tolerance:69,(+)Tolerance:71 Standard Weight:80,(-)Tolerance:79,(+)Tolerance:81	d

AUDIT TRAIL FOR BALANCE

DATE	TIME	USER ID	ACTION	BALANCE CODE NO.	OLD VALUES	NEW VALUES	OLD WEIGHT	NEW WEIGHT	REMARK
13.10.2024	10:13:37	User ID:123,User Name:yogita	Add	Balance-121212	NA	Balance Code No.: NA Balance-121212, Make: METTLER TOLEDO, Model: ML303, Serial No.: NA, Unit: gm, Least Count: 0.001, Maximum Capacity: 350.000, Minimum Capacity: 0.001, Maximum Operating Range: 250.000, Minimum Operating Range: 0.002, Periodic Data: Calibration Date:13-10-2024, Duration :5, Reminder:4, Balance Type:Analytical	NA	Daily: Standard Weight:0.020,(-)Tolerance:0.018,(+)Tolerance:0.022 Standard Weight:10.000,(-)Tolerance:9.896,(+)Tolerance:10.452 Standard Weight:20.000,(-)Tolerance:18.000,(+)Tolerance:22.000 Standard Weight:30.000,(-)Tolerance:28.000,(+)Tolerance:31.000 0 Eccentricity: Standard Weight:0.020,(-)Tolerance:0.018,(+)Tolerance:0.022 Repeatability: Standard Weight:0.020,(-)Tolerance:0.018,(+)Tolerance:0.022 Periodic: Standard Weight:0.020,(-)Tolerance:0.018,(+)Tolerance:0.022 Standard Weight:10.000,(-)Tolerance:9.896,(+)Tolerance:10.452 Standard Weight:20.000,(-)Tolerance:18.000,(+)Tolerance:22.000 Standard Weight:30.000,(-)Tolerance:28.000,(+)Tolerance:31.000 Standard Weight:40.000,(-)Tolerance:39.000,(+)Tolerance:41.000 Standard Weight:50.000,(-)Tolerance:48.000,(+)Tolerance:51.000 Standard Weight:60.000,(-)Tolerance:58.000,(+)Tolerance:61.000 0 Sensitivity: Standard Weight:0.020,(-)Tolerance:0.018,(+)Tolerance:0.022	6