

AUDIT TRAIL FOR BALANCE

DATE	TIME	USER ID	ACTION	BALANCE ID	OLD VALUES	NEW VALUES	OLD WEIGHT	NEW WEIGHT	REMARK
18.09.2024	14:23:56	User ID:1,User Name:Ashwin Reddy	Activate	Balance3	NA	NA	NA	NA	as

AUDIT TRAIL FOR BALANCE

DATE	TIME	USER ID	ACTION	BALANCE ID	OLD VALUES	NEW VALUES	OLD WEIGHT	NEW WEIGHT	REMARK
18.09.2024	14:22:44	User ID:11,User Name:lavya singh	Add	Balance3	NA	Balance ID: Balance3,Make: MT,Model: PB303S,Serial No.: Serial3,Balance Certificate No.:NA,Unit: gm,Least Count: 0.1,Maximum Capacity: 100,Minimum Capacity: 0.2,Maximum Operating Range: 99,Minimum Operating Range: 0.3,Periodic Data: Calibration Date:16-10- 2024,Duration :5,Reminder:2,Balance Type:Analytical	NA	Daily: Standard Weight:5,(-)Tolerance:4,(+)Tolerance:6 Standard Weight:10,(-)Tolerance:9,(+)Tolerance:11 Standard Weight:15,(-)Tolerance:14,(+)Tolerance:16 Standard Weight:20,(-)Tolerance:19,(+)Tolerance:21 Periodic: Standard Weight:5,(-)Tolerance:4,(+)Tolerance:6 Standard Weight:10,(-)Tolerance:9,(+)Tolerance:11 Standard Weight:15,(-)Tolerance:14,(+)Tolerance:16 Standard Weight:20,(-)Tolerance:19,(+)Tolerance:21 Standard Weight:25,(-)Tolerance:24,(+)Tolerance:26 Standard Weight:30,(-)Tolerance:29,(+)Tolerance:31 Standard Weight:35,(-)Tolerance:34,(+)Tolerance:36 Standard Weight:40,(-)Tolerance:39,(+)Tolerance:41 Standard Weight:45,(-)Tolerance:44,(+)Tolerance:46 Standard Weight:50,(-)Tolerance:49,(+)Tolerance:51	as

Page 3 of 6

AUDIT TRAIL FOR BALANCE

DATE	TIME	USER ID	ACTION	BALANCE ID	OLD VALUES	NEW VALUES	OLD WEIGHT	NEW WEIGHT	REMARK
18.09.2024	12:36:26	User ID:1,User Name:Ashwin Reddy	Add	IPQCBalance1	NA	Balance ID: IPQCBalance1,Make: MT,Model: PB303S,Serial No.: IPQCSeria1,Balance Certificate No.:NA,Unit: gm,Least Count: 0.1,Maximum Capacity: 100,Minimum Capacity: 0.2,Maximum Operating Range: 99,Minimum Operating Range: 0.3,Periodic Data: Calibration Date:27-09- 2024,Duration :5,Reminder:2,Balance Type:Analytical	NA	Daily: Standard Weight:5,(-)Tolerance:4,(+)Tolerance:6 Standard Weight:10,(-)Tolerance:9,(+)Tolerance:11 Standard Weight:15,(-)Tolerance:14,(+)Tolerance:16 Standard Weight:20,(-)Tolerance:19,(+)Tolerance:21 Periodic: Standard Weight:5,(-)Tolerance:4,(+)Tolerance:6 Standard Weight:10,(-)Tolerance:9,(+)Tolerance:11 Standard Weight:15,(-)Tolerance:14,(+)Tolerance:16 Standard Weight:20,(-)Tolerance:19,(+)Tolerance:21 Standard Weight:25,(-)Tolerance:24,(+)Tolerance:26 Standard Weight:30,(-)Tolerance:29,(+)Tolerance:31 Standard Weight:35,(-)Tolerance:34,(+)Tolerance:36 Standard Weight:40,(-)Tolerance:39,(+)Tolerance:41 Standard Weight:45,(-)Tolerance:44,(+)Tolerance:46 Standard Weight:50,(-)Tolerance:49,(+)Tolerance:51	as

AUDIT TRAIL FOR BALANCE

DATE	TIME	USER ID	ACTION	BALANCE ID	OLD VALUES	NEW VALUES	OLD WEIGHT	NEW WEIGHT	REMARK
18.09.2024	10:01:25	User ID:11,User Name:lavya singh	Activate	Balance1	NA	NA	NA	NA	as

AUDIT TRAIL FOR BALANCE

DATE	TIME	USER ID	ACTION	BALANCE ID	OLD VALUES	NEW VALUES	OLD WEIGHT	NEW WEIGHT	REMARK
18.09.2024	09:59:04	User ID:1,User Name:Ashwin Reddy	Add	Balance1	NA	Balance ID: Balance1,Make: MT,Model: PB303S,Serial No.: Serial1,Balance Certificate No.:NA,Unit: gm,Least Count: 0.1,Maximum Capacity: 100,Minimum Capacity: 0.2,Maximum Operating Range: 99,Minimum Operating Range: 0.3,Periodic Data: Calibration Date:10-10- 2024,Duration :5,Reminder:2,Balance Type:Analytical	NA	Daily: Standard Weight:5,(-)Tolerance:4,(+)Tolerance:6 Standard Weight:10,(-)Tolerance:9,(+)Tolerance:11 Standard Weight:15,(-)Tolerance:14,(+)Tolerance:16 Standard Weight:20,(-)Tolerance:19,(+)Tolerance:21 Periodic: Standard Weight:5,(-)Tolerance:4,(+)Tolerance:6 Standard Weight:10,(-)Tolerance:9,(+)Tolerance:11 Standard Weight:15,(-)Tolerance:14,(+)Tolerance:16 Standard Weight:20,(-)Tolerance:19,(+)Tolerance:21 Standard Weight:25,(-)Tolerance:24,(+)Tolerance:26 Standard Weight:30,(-)Tolerance:29,(+)Tolerance:31 Standard Weight:35,(-)Tolerance:34,(+)Tolerance:36 Standard Weight:40,(-)Tolerance:39,(+)Tolerance:41 Standard Weight:45,(-)Tolerance:44,(+)Tolerance:46 Standard Weight:50,(-)Tolerance:49,(+)Tolerance:51	as