

FULL SCALE CALIBRATION AND ZERO ERROR RECORD
(Analytical and Precision Weighing Balance)

Unit : X
Department : QC
Balance Code No. : BCN/3/BAL/09
Make : Sartorius
Std. Weight Box ID. : WB_10

Month and year : January 2026
Capacity (a) (g) : 356.000
Model : BC-ED
Least Count (e) (g) 0.001

Std. Weight Box Certificate No. 240604

Std. Weight Box Valid Upto 17-08-2030

Tolerance: +/- Least count of the balance OR +/- 0.1 % of Standard weight whichever is higher.

Standard weights for calibration (g)		Tolerance (g)	Accuracy Limits(Std wt +/- Tolerance) (g)
Lower(20.e)/(5% of a)/Minimum capacity of the checkweigher / weighing balance	1.970	- 0.470 + 0.380	1.500 - 2.350
Middle-I(20% of a)	4.970	- 1.470 + 0.280	3.500 - 5.250
Middle-II(50% of a)	9.920	- 0.470 + 0.230	9.450 - 10.150
Upper(80% of a)	19.940	- 0.710 + 0.560	19.230 - 20.500

Note:

1. Adjust zero error before starting the weighing operation.
2. Check against standard weights daily before starting the activity.
3. If more than one standard weight is used for verification then record the standard weight / identification No. of all individual weights used for verification. (E.g. For verification of 200g, if three weights viz. 50g, 50g & 100g is used then record standard weight/identification No. of all the three weights)

Date	Time	Spirit Level (Ok/ Not Ok)	Zero Error	Standard weights kept on the balance (g)		Readings shown (g)	Remark
13.01.2026	11:42	OK	Nil	Lower(20.e)/(5% of a)/Minimum capacity of the checkweigher / weighing balance	1.970	1.970	Ok
				Middle-I(20% of a)	4.970	5.200	Ok
				Middle-II(50% of a)	9.920	9.952	Ok
				Upper(80% of a)	19.940	19.945	Ok

	Signature	Date & Time
Performed By	TANTRA (TANTRA)	13.01.2026 & 11:42:05
Checked By		
Checked By (Quality Assurance)		